



The Clipper. **Redesigned.**
Power. Comfort. Safety.

The new Medline clipper

Design with purpose

High grade steel blades | Powerful motor | Extensive usability testing

Designed for Performance

Head-to-head testing by an independent lab confirmed Medline's universal blade exhibited equivalent performance against the industry leader¹

Designed for Patients

- Single-patient-use disposable blades for sharpness and safety
- Moving blade kept away from patient's skin to reduce nicks and cuts
- Used wet or dry, effective on multiple hair types

Designed for Simplicity

- Simple, one-switch control
- Wall-mount or freestanding charging base

Designed for Outcomes

- The clip approach favored by AORN^{2,3} and SCIP⁴
- Water-resistant—easy to disinfect

Designed for Clinicians

- "Grip-like-a-pen" ergonomics for comfort and control
- Easy-on / Easy-off single-patient use disposable blades



Designed for a natural grip and optimal cutting angle



Easy-on/easy-off clipper head with alignment marks for secure attachment



Designed for Versatility

Disposable blades available for sensitive areas, general purpose and coarse hair types



Works wet or dry for simple cleaning, disinfection and use



Designed for Value

The switch to the new Medline clipper can be more economical than you'd imagine, whether you're using a competitive clipper or traditional razor. Your Medline rep has all the attractive details about our clipper program.



Item #	Description	Pkg.
DYND70800	Medline Surgical Clipper, Charger and Base	1 ea/bx
DYND70801	Medline Surgical Clipper, Handle Only	1 ea/bx
DYND70880	Universal Body Hair Clipping Blade, Gray	50 ea/cs
DYND70885	Coarse Hair Clipping Blade, Blue	50 ea/cs
DYND70890	Narrow/Sensitive Area Clipping Blade, Red	50 ea/cs

DYND70800

Charger stands alone or mounts to wall



DYND70890



DYND70880



DYND70885

References

1. LNE test report, June 2013. Data on file, available on request.
2. Perioperative Standards and Recommended Practices. Denver, CO: AORN, Inc; 2013:75-90. Accessed at <http://www.aorn.org/Secondary.aspx?id=20981> November 2013.
3. Kjonniksen I, Andersen BM, Sondenaa VG, et al. Preoperative hair removal-a systematic literature review. AORN J. 2002 May;75 (5):928-938,940. PMID:12063942.
4. Specifications Manual for National Hospital Inpatient Quality Measures (SCIP-Inf 6) Discharges 10-01-09 (4Q09) through 03-31-10 (1Q10).



Medline Industries, Inc.
One Medline Place
Mundelein, IL 60060

Medline United States
1-800-MEDLINE (633-5463)
www.medline.com
info@medline.com

Medline Canada
1-800-396-6996
www.medline.ca
canada@medline.com

Medline México
01-800-831-0898
www.medlinemexico.com
mexico@medline.com

Some products may not be available for sale in Mexico or Canada. We reserve the right to correct any errors that may occur within this brochure. ©2013 Medline Industries, Inc. All rights reserved. Medline is a registered trademark of Medline Industries, Inc.