

SUCTION THERAPY

Suction Regulators

Solutions for Life

88XVR Series Suction Regulators

Mode Selector Switch

- Clearly marked indicator
- GENTEC® suction regulators are available in continuous only (882 Series) and continuous-intermittent (881 Series)
- 160 mmHg and 200 mmHg are available in pediatric, 300 mmHg is available in adult, and 760 mmHg is available in surgical
- Factory default settings for INT (on / off) cycles on the 881 Series is 8 seconds off, 16 seconds on. On/off settings can be adjusted

Control Knob

- Easy adjustment of vacuum pressure
- No Leakage
- Agion® Antimicrobial protection
 - Built-in long lasting antimicrobial protection
 - Based on naturally-occurring silver
 - Inhibits growth of microbes on surface of knobs and switches
 - Environmentally friendly technology
 - Broad Regulatory Approvals including EPA, FDA, NSF, USDA, Canada, European BPD, Swiss and EFSA
- Color Coding available in both USA and ISO standard



Fittings

- GENTEC® provides a wide range of fittings to meet customers' needs



881VR-300D-XX
Continuous / Intermittent
0-300 mmHg, digital gauge

Glow in the Dark Gauge

- 2-1/2" gauge, easy to read
- Color coded flow ranges
- Dual Scale mmHg and kPa
- Clockwise, counter-clockwise or digital gauge mode (features an analog scale in digital mode where a digital reading is on screen and an analog meter on the outer ring)



Face Cover

- High-quality ABS case body for long durability
- Impact-resistant internal components minimize damage if regulator is bumped or dropped
- Easy to clean, wipe down with sterilant

Interior Structure


















- Few mechanical parts for easy maintenance
- Corrosion free
- Lightweight

Suction Trap

- Prevents contamination caused by overflow



* Agion's silver antimicrobial has been certified for its material content, recyclability, and manufacturing characteristics. Cradle To Cradle Certified CM products meet established standards for human health, environmental health and recyclability. Cradle To Cradle Certified CM is a certification mark of MBDC.

	881VR Series Continuous / Intermittent			882VR Series Continuous-only			
Clockwise	 881VR-160 pediatric	 881VR-200 adult	 881VR-300 adult	 882VR-160 pediatric	 882VR-200 adult	 882VR-300 adult	 882VR-760 surgical
Counter-Clockwise	 881VR-160E pediatric		 881VR-300E adult	 882VR-160E adult		 882VR-300E adult	 882VR-760E surgical
Digital Gauge	 881VR-160D pediatric		 881VR-300D adult	 882VR-160D pediatric		 882VR-300D adult	 882VR-760D surgical
Mode Selection:							
OFF	√	√	√	√	√	√	√
CONT	√	√	√	√	√	√	√
INT	√	√	√	X	X	X	X
FULL	X	X	X	X	√	√	√
Regulated Vacuum Range	0~160 mmHg	0~200 mmHg	0~300 mmHg	0~160 mmHg	0~200 mmHg	0~300 mmHg	0~760 mmHg
Gauge Display Range	0~160 mmHg	0~200 mmHg	0~300 mmHg	0~160 mmHg	0~200 mmHg	0~300 mmHg	0~760 mmHg
Gauge Accuracy	± 4%	± 4%	± 4%	± 4%	± 4%	± 4%	± 4%
Intermittent Time	Suction: 16 seconds, Stop: 8 seconds			/	/	/	/
Flow Rate	>50 l/min (CONT) 7-9 l/min (INT)	>55 l/min (CONT) 7-9 l/min (INT)	>60 l/min (CONT) 7-9 l/min (INT)	>50 l/min	>55 l/min	>60 l/min	>70 l/min
Weight	1 lb / 0.455 kg			0.81 lb / 0.37 kg			
Dimensions	6.3" H × 3.15" W × 3.94" D / 16 × 8 × 10 cm			6.3" H × 3.15" W × 3.94" D / 16 × 8 × 10 cm			

SUCTION THERAPY

Suction Regulators

Solutions for Life

88XVR Series Suction Regulators (Range 160)

 <p>Vacuum Connector</p> <p>Suction Connector</p>	 <p>No Adapter 1/8" NPT (F)</p>	 <p>OH: Omeda® 3111-VAC-NT2</p>	 <p>DH: DISS Hand-tight 3222-VAC-NT2</p>	 <p>DS: DISS Hex 3102-VAC-NT2</p>	 <p>CH: Chemetron® 3103-VAC-NT2</p>
 <p>No Adapter 1/8" NPT (F)</p>	<p>881VR-160-00-A 881VR-160E-00-A 881VR-160D-00-A 882VR-160-00-A 882VR-160E-00-A 882VR-160D-00-A 882VRM-160-00-A 882VRM-160E-00-A</p>	<p>881VR-160-OH-A 881VR-160E-OH-A 881VR-160D-OH-A 882VR-160-OH-A 882VR-160E-OH-A 882VR-160D-OH-A 882VRM-160-OH-A 882VRM-160E-OH-A</p>	<p>881VR-160-DH-A 881VR-160E-DH-A 881VR-160D-DH-A 882VR-160-DH-A 882VR-160E-DH-A 882VR-160D-DH-A 882VRM-160-DH-A 882VRM-160E-DH-A</p>	<p>881VR-160-DS-A 881VR-160E-DS-A 881VR-160D-DS-A 882VR-160-DS-A 882VR-160E-DS-A 882VR-160D-DS-A 882VRM-160-DS-A 882VRM-160E-DS-A</p>	<p>881VR-160-CH-A 881VR-160E-CH-A 881VR-160D-CH-A 882VR-160-CH-A 882VR-160E-CH-A 882VR-160D-CH-A 882VRM-160-CH-A 882VRM-160E-CH-A</p>
 <p>Suction Trap Assembly (complete)</p>	<p>881VR-160-00-T 881VR-160E-00-T 881VR-160D-00-T 882VR-160-00-T 882VR-160E-00-T 882VR-160D-00-T 882VRM-160-00-T 882VRM-160E-00-T</p>	<p>881VR-160-OH-T 881VR-160E-OH-T 881VR-160D-OH-T 882VR-160-OH-T 882VR-160E-OH-T 882VR-160D-OH-T 882VRM-160-OH-T 882VRM-160E-OH-T</p>	<p>881VR-160-DH-T 881VR-160E-DH-T 881VR-160D-DH-T 882VR-160-DH-T 882VR-160E-DH-T 882VR-160D-DH-T 882VRM-160-DH-T 882VRM-160E-DH-T</p>	<p>881VR-160-DS-T 881VR-160E-DS-T 881VR-160D-DS-T 882VR-160-DS-T 882VR-160E-DS-T 882VR-160D-DS-T 882VRM-160-DS-T 882VRM-160E-DS-T</p>	<p>881VR-160-CH-T 881VR-160E-CH-T 881VR-160D-CH-T 882VR-160-CH-T 882VR-160E-CH-T 882VR-160D-CH-T 882VRM-160-CH-T 882VRM-160E-CH-T</p>
 <p>Trap Fitting Only (No Trap)</p>	<p>881VR-160-00-B 881VR-160E-00-B 881VR-160D-00-B 882VR-160-00-B 882VR-160E-00-B 882VR-160D-00-B 882VRM-160-00-B 882VRM-160E-00-B</p>	<p>881VR-160-OH-B 881VR-160E-OH-B 881VR-160D-OH-B 882VR-160-OH-B 882VR-160E-OH-B 882VR-160D-OH-B 882VRM-160-OH-B 882VRM-160E-OH-B</p>	<p>881VR-160-DH-B 881VR-160E-DH-B 881VR-160D-DH-B 882VR-160-DH-B 882VR-160E-DH-B 882VR-160D-DH-B 882VRM-160-DH-B 882VRM-160E-DH-B</p>	<p>881VR-160-DS-B 881VR-160E-DS-B 881VR-160D-DS-B 882VR-160-DS-B 882VR-160E-DS-B 882VR-160D-DS-B 882VRM-160-DS-B 882VRM-160E-DS-B</p>	<p>881VR-160-CH-B 881VR-160E-CH-B 881VR-160D-CH-B 882VR-160-CH-B 882VR-160E-CH-B 882VR-160D-CH-B 882VRM-160-CH-B 882VRM-160E-CH-B</p>
 <p>DISS Male Vacuum</p>	<p>881VR-160-00-C 881VR-160E-00-C 881VR-160D-00-C 882VR-160-00-C 882VR-160E-00-C 882VR-160D-00-C 882VRM-160-00-C 882VRM-160E-00-C</p>	<p>881VR-160-OH-C 881VR-160E-OH-C 881VR-160D-OH-C 882VR-160-OH-C 882VR-160E-OH-C 882VR-160D-OH-C 882VRM-160-OH-C 882VRM-160E-OH-C</p>	<p>881VR-160-DH-C 881VR-160E-DH-C 881VR-160D-DH-C 882VR-160-DH-C 882VR-160E-DH-C 882VR-160D-DH-C 882VRM-160-DH-C 882VRM-160E-DH-C</p>	<p>881VR-160-DS-C 881VR-160E-DS-C 881VR-160D-DS-C 882VR-160-DS-C 882VR-160E-DS-C 882VR-160D-DS-C 882VRM-160-DS-C 882VRM-160E-DS-C</p>	<p>881VR-160-CH-C 881VR-160E-CH-C 881VR-160D-CH-C 882VR-160-CH-C 882VR-160E-CH-C 882VR-160D-CH-C 882VRM-160-CH-C 882VRM-160E-CH-C</p>
 <p>Tubing Nipple (5/16" Hose Barb)</p>	<p>881VR-160-00-D 881VR-160E-00-D 881VR-160D-00-D 882VR-160-00-D 882VR-160E-00-D 882VR-160D-00-D 882VRM-160-00-D 882VRM-160E-00-D</p>	<p>881VR-160-OH-D 881VR-160E-OH-D 881VR-160D-OH-D 882VR-160-OH-D 882VR-160E-OH-D 882VR-160D-OH-D 882VRM-160-OH-D 882VRM-160E-OH-D</p>	<p>881VR-160-DH-D 881VR-160E-DH-D 881VR-160D-DH-D 882VR-160-DH-D 882VR-160E-DH-D 882VR-160D-DH-D 882VRM-160-DH-D 882VRM-160E-DH-D</p>	<p>881VR-160-DS-D 881VR-160E-DS-D 881VR-160D-DS-D 882VR-160-DS-D 882VR-160E-DS-D 882VR-160D-DS-D 882VRM-160-DS-D 882VRM-160E-DS-D</p>	<p>881VR-160-CH-D 881VR-160E-CH-D 881VR-160D-CH-D 882VR-160-CH-D 882VR-160E-CH-D 882VR-160D-CH-D 882VRM-160-CH-D 882VRM-160E-CH-D</p>

 <p>PB: Puritan-Bennett® 3114-VAC-NT2</p>	 <p>FS 3101-VAC-NT2</p>	 <p>GS 3116-VAC-NT2</p>	 <p>BS 3117-VAC-NT2</p>	 <p>JIS 3108-VAC-NT2</p>	 <p>SIS 3109-VAC-NT2</p>
881VR-160-PB-A 881VR-160E-PB-A 881VR-160D-PB-A 882VR-160-PB-A 882VR-160E-PB-A 882VR-160D-PB-A 882VRM-160-PB-A 882VRM-160E-PB-A	881VR-160-FS-A 881VR-160E-FS-A 881VR-160D-FS-A 882VR-160-FS-A 882VR-160E-FS-A 882VR-160D-FS-A 882VRM-160-FS-A 882VRM-160E-FS-A	881VR-160-GS-A 881VR-160E-GS-A 881VR-160D-GS-A 882VR-160-GS-A 882VR-160E-GS-A 882VR-160D-GS-A 882VRM-160-GS-A 882VRM-160E-GS-A	881VR-160-BS-A 881VR-160E-BS-A 881VR-160D-BS-A 882VR-160-BS-A 882VR-160E-BS-A 882VR-160D-BS-A 882VRM-160-BS-A 882VRM-160E-BS-A	881VR-160-JIS-A 881VR-160E-JIS-A 881VR-160D-JIS-A 882VR-160-JIS-A 882VR-160E-JIS-A 882VR-160D-JIS-A 882VRM-160-JIS-A 882VRM-160E-JIS-A	881VR-160-SIS-A 881VR-160E-SIS-A 881VR-160D-SIS-A 882VR-160-SIS-A 882VR-160E-SIS-A 882VR-160D-SIS-A 882VRM-160-SIS-A 882VRM-160E-SIS-A
881VR-160-PB-T 881VR-160E-PB-T 881VR-160D-PB-T 882VR-160-PB-T 882VR-160E-PB-T 882VR-160D-PB-T 882VRM-160-PB-T 882VRM-160E-PB-T	881VR-160-FS-T 881VR-160E-FS-T 881VR-160D-FS-T 882VR-160-FS-T 882VR-160E-FS-T 882VR-160D-FS-T 882VRM-160-FS-T 882VRM-160E-FS-T	881VR-160-GS-T 881VR-160E-GS-T 881VR-160D-GS-T 882VR-160-GS-T 882VR-160E-GS-T 882VR-160D-GS-T 882VRM-160-GS-T 882VRM-160E-GS-T	881VR-160-BS-T 881VR-160E-BS-T 881VR-160D-BS-T 882VR-160-BS-T 882VR-160E-BS-T 882VR-160D-BS-T 882VRM-160-BS-T 882VRM-160E-BS-T	881VR-160-JIS-T 881VR-160E-JIS-T 881VR-160D-JIS-T 882VR-160-JIS-T 882VR-160E-JIS-T 882VR-160D-JIS-T 882VRM-160-JIS-T 882VRM-160E-JIS-T	881VR-160-SIS-T 881VR-160E-SIS-T 881VR-160D-SIS-T 882VR-160-SIS-T 882VR-160E-SIS-T 882VR-160D-SIS-T 882VRM-160-SIS-T 882VRM-160E-SIS-T
881VR-160-PB-B 881VR-160E-PB-B 881VR-160D-PB-B 882VR-160-PB-B 882VR-160E-PB-B 882VR-160D-PB-B 882VRM-160-PB-B 882VRM-160E-PB-B	881VR-160-FS-B 881VR-160E-FS-B 881VR-160D-FS-B 882VR-160-FS-B 882VR-160E-FS-B 882VR-160D-FS-B 882VRM-160-FS-B 882VRM-160E-FS-B	881VR-160-GS-B 881VR-160E-GS-B 881VR-160D-GS-B 882VR-160-GS-B 882VR-160E-GS-B 882VR-160D-GS-B 882VRM-160-GS-B 882VRM-160E-GS-B	881VR-160-BS-B 881VR-160E-BS-B 881VR-160D-BS-B 882VR-160-BS-B 882VR-160E-BS-B 882VR-160D-BS-B 882VRM-160-BS-B 882VRM-160E-BS-B	881VR-160-JIS-B 881VR-160E-JIS-B 881VR-160D-JIS-B 882VR-160-JIS-B 882VR-160E-JIS-B 882VR-160D-JIS-B 882VRM-160-JIS-B 882VRM-160E-JIS-B	881VR-160-SIS-B 881VR-160E-SIS-B 881VR-160D-SIS-B 882VR-160-SIS-B 882VR-160E-SIS-B 882VR-160D-SIS-B 882VRM-160-SIS-B 882VRM-160E-SIS-B
881VR-160-PB-C 881VR-160E-PB-C 881VR-160D-PB-C 882VR-160-PB-C 882VR-160E-PB-C 882VR-160D-PB-C 882VRM-160-PB-C 882VRM-160E-PB-C	881VR-160-FS-C 881VR-160E-FS-C 881VR-160D-FS-C 882VR-160-FS-C 882VR-160E-FS-C 882VR-160D-FS-C 882VRM-160-FS-C 882VRM-160E-FS-C	881VR-160-GS-C 881VR-160E-GS-C 881VR-160D-GS-C 882VR-160-GS-C 882VR-160E-GS-C 882VR-160D-GS-C 882VRM-160-GS-C 882VRM-160E-GS-C	881VR-160-BS-C 881VR-160E-BS-C 881VR-160D-BS-C 882VR-160-BS-C 882VR-160E-BS-C 882VR-160D-BS-C 882VRM-160-BS-C 882VRM-160E-BS-C	881VR-160-JIS-C 881VR-160E-JIS-C 881VR-160D-JIS-C 882VR-160-JIS-C 882VR-160E-JIS-C 882VR-160D-JIS-C 882VRM-160-JIS-C 882VRM-160E-JIS-C	881VR-160-SIS-C 881VR-160E-SIS-C 881VR-160D-SIS-C 882VR-160-SIS-C 882VR-160E-SIS-C 882VR-160D-SIS-C 882VRM-160-SIS-C 882VRM-160E-SIS-C
881VR-160-PB-D 881VR-160E-PB-D 881VR-160D-PB-D 882VR-160-PB-D 882VR-160E-PB-D 882VR-160D-PB-D 882VRM-160-PB-D 882VRM-160E-PB-D	881VR-160-FS-D 881VR-160E-FS-D 881VR-160D-FS-D 882VR-160-FS-D 882VR-160E-FS-D 882VR-160D-FS-D 882VRM-160-FS-D 882VRM-160E-FS-D	881VR-160-GS-D 881VR-160E-GS-D 881VR-160D-GS-D 882VR-160-GS-D 882VR-160E-GS-D 882VR-160D-GS-D 882VRM-160-GS-D 882VRM-160E-GS-D	881VR-160-BS-D 881VR-160E-BS-D 881VR-160D-BS-D 882VR-160-BS-D 882VR-160E-BS-D 882VR-160D-BS-D 882VRM-160-BS-D 882VRM-160E-BS-D	881VR-160-JIS-D 881VR-160E-JIS-D 881VR-160D-JIS-D 882VR-160-JIS-D 882VR-160E-JIS-D 882VR-160D-JIS-D 882VRM-160-JIS-D 882VRM-160E-JIS-D	881VR-160-SIS-D 881VR-160E-SIS-D 881VR-160D-SIS-D 882VR-160-SIS-D 882VR-160E-SIS-D 882VR-160D-SIS-D 882VRM-160-SIS-D 882VRM-160E-SIS-D

see
previous
page for
complete
table

SUCTION THERAPY

Suction Regulators

Solutions for Life

88XVR Series Suction Regulators (Range 200)

 <p>Vacuum Connector</p> <p>Suction Connector</p>	 <p>No Adapter 1/8" NPT (F)</p>	 <p>OH: Omeda® 3111-VAC-NT2</p>	 <p>DH: DISS Hand-tight 3222-VAC-NT2</p>	 <p>DS: DISS Hex 3102-VAC-NT2</p>	 <p>CH: Chemetron® 3103-VAC-NT2</p>
 <p>No Adapter 1/8" NPT (F)</p>	<p>881VR-200-00-A 882VR-200-00-A 882VRM-200-00-A</p>	<p>881VR-200-OH-A 882VR-200-OH-A 882VRM-200-OH-A</p>	<p>881VR-200-DH-A 882VR-200-DH-A 882VRM-200-DH-A</p>	<p>881VR-200-DS-A 882VR-200-DS-A 882VRM-200-DS-A</p>	<p>881VR-200-CH-A 882VR-200-CH-A 882VRM-200-CH-A</p>
 <p>Suction Trap Assembly (complete)</p>	<p>881VR-200-00-T 882VR-200-00-T 882VRM-200-00-T</p>	<p>881VR-200-OH-T 882VR-200-OH-T 882VRM-200-OH-T</p>	<p>881VR-200-DH-T 882VR-200-DH-T 882VRM-200-DH-T</p>	<p>881VR-200-DS-T 882VR-200-DS-T 882VRM-200-DS-T</p>	<p>881VR-200-CH-T 882VR-200-CH-T 882VRM-200-CH-T</p>
 <p>Trap Fitting Only (No Trap)</p>	<p>881VR-200-00-B 882VR-200-00-B 882VRM-200-00-B</p>	<p>881VR-200-OH-B 882VR-200-OH-B 882VRM-200-OH-B</p>	<p>881VR-200-DH-B 882VR-200-DH-B 882VRM-200-DH-B</p>	<p>881VR-200-DS-B 882VR-200-DS-B 882VRM-200-DS-B</p>	<p>881VR-200-CH-B 882VR-200-CH-B 882VRM-200-CH-B</p>
 <p>DISS Male Vacuum</p>	<p>881VR-200-00-C 882VR-200-00-C 882VRM-200-00-C</p>	<p>881VR-200-OH-C 882VR-200-OH-C 882VRM-200-OH-C</p>	<p>881VR-200-DH-C 882VR-200-DH-C 882VRM-200-DH-C</p>	<p>881VR-200-DS-C 882VR-200-DS-C 882VRM-200-DS-C</p>	<p>881VR-200-CH-C 882VR-200-CH-C 882VRM-200-CH-C</p>
 <p>Tubing Nipple (5/16" Hose Barb)</p>	<p>881VR-200-00-D 882VR-200-00-D 882VRM-200-00-D</p>	<p>881VR-200-OH-D 882VR-200-OH-D 882VRM-200-OH-D</p>	<p>881VR-200-DH-D 882VR-200-DH-D 882VRM-200-DH-D</p>	<p>881VR-200-DS-D 882VR-200-DS-D 882VRM-200-DS-D</p>	<p>881VR-200-CH-D 882VR-200-CH-D 882VRM-200-CH-D</p>

 <p>PB: Puritan-Bennett® 3114-VAC-NT2</p>	 <p>FS 3101-VAC-NT2</p>	 <p>GS 3116-VAC-NT2</p>	 <p>BS 3117-VAC-NT2</p>	 <p>JIS 3108-VAC-NT2</p>	 <p>SIS 3109-VAC-NT2</p>
<p>881VR-200-PB-A 882VR-200-PB-A 882VRM-200-PB-A</p>	<p>881VR-200-FS-A 882VR-200-FS-A 882VRM-200-FS-A</p>	<p>881VR-200-GS-A 882VR-200-GS-A 882VRM-200-GS-A</p>	<p>881VR-200-BS-A 882VR-200-BS-A 882VRM-200-BS-A</p>	<p>881VR-200-JIS-A 882VR-200-JIS-A 882VRM-200-JIS-A</p>	<p>881VR-200-SIS-A 882VR-200-SIS-A 882VRM-200-SIS-A</p>
<p>881VR-200-PB-T 882VR-200-PB-T 882VRM-200-PB-T</p>	<p>881VR-200-FS-T 882VR-200-FS-T 882VRM-200-FS-T</p>	<p>881VR-200-GS-T 882VR-200-GS-T 882VRM-200-GS-T</p>	<p>881VR-200-BS-T 882VR-200-BS-T 882VRM-200-BS-T</p>	<p>881VR-200-JIS-T 882VR-200-JIS-T 882VRM-200-JIS-T</p>	<p>881VR-200-SIS-T 882VR-200-SIS-T 882VRM-200-SIS-T</p>
<p>881VR-200-PB-B 882VR-200-PB-B 882VRM-200-PB-B</p>	<p>881VR-200-FS-B 882VR-200-FS-B 882VRM-200-FS-B</p>	<p>881VR-200-GS-B 882VR-200-GS-B 882VRM-200-GS-B</p>	<p>881VR-200-BS-B 882VR-200-BS-B 882VRM-200-BS-B</p>	<p>881VR-200-JIS-B 882VR-200-JIS-B 882VRM-200-JIS-B</p>	<p>881VR-200-SIS-B 882VR-200-SIS-B 882VRM-200-SIS-B</p>
<p>881VR-200-PB-C 882VR-200-PB-C 882VRM-200-PB-C</p>	<p>881VR-200-FS-C 882VR-200-FS-C 882VRM-200-FS-C</p>	<p>881VR-200-GS-C 882VR-200-GS-C 882VRM-200-GS-C</p>	<p>881VR-200-BS-C 882VR-200-BS-C 882VRM-200-BS-C</p>	<p>881VR-200-JIS-C 882VR-200-JIS-C 882VRM-200-JIS-C</p>	<p>881VR-200-SIS-C 882VR-200-SIS-C 882VRM-200-SIS-C</p>
<p>881VR-200-PB-D 882VR-200-PB-D 882VRM-200-PB-D</p>	<p>881VR-200-FS-D 882VR-200-FS-D 882VRM-200-FS-D</p>	<p>881VR-200-GS-D 882VR-200-GS-D 882VRM-200-GS-D</p>	<p>881VR-200-BS-D 882VR-200-BS-D 882VRM-200-BS-D</p>	<p>881VR-200-JIS-D 882VR-200-JIS-D 882VRM-200-JIS-D</p>	<p>881VR-200-SIS-D 882VR-200-SIS-D 882VRM-200-SIS-D</p>

see
previous
page for
complete
table

SUCTION THERAPY

Suction Regulators

Solutions for Life

88XVR Series Suction Regulators (Range 300)

 <p>Vacuum Connector</p> <p>Suction Connector</p>	 <p>No Adapter 1/8" NPT (F)</p>	 <p>OH: Omeda® 3111-VAC-NT2</p>	 <p>DH: DISS Hand-tight 3222-VAC-NT2</p>	 <p>DS: DISS Hex 3102-VAC-NT2</p>	 <p>CH: Chemetron® 3103-VAC-NT2</p>
 <p>No Adapter 1/8" NPT (F)</p>	<p>881VR-300-00-A 881VR-300E-00-A 881VR-300D-00-A 882VR-300-00-A 882VR-300E-00-A 882VR-300D-00-A 882VRM-300-00-A 882VRM-300E-00-A</p>	<p>881VR-300-OH-A 881VR-300E-OH-A 881VR-300D-OH-A 882VR-300-OH-A 882VR-300E-OH-A 882VR-300D-OH-A 882VRM-300-OH-A 882VRM-300E-OH-A</p>	<p>881VR-300-DH-A 881VR-300E-DH-A 881VR-300D-DH-A 882VR-300-DH-A 882VR-300E-DH-A 882VR-300D-DH-A 882VRM-300-DH-A 882VRM-300E-DH-A</p>	<p>881VR-300-DS-A 881VR-300E-DS-A 881VR-300D-DS-A 882VR-300-DS-A 882VR-300E-DS-A 882VR-300D-DS-A 882VRM-300-DS-A 882VRM-300E-DS-A</p>	<p>881VR-300-CH-A 881VR-300E-CH-A 881VR-300D-CH-A 882VR-300-CH-A 882VR-300E-CH-A 882VR-300D-CH-A 882VRM-300-CH-A 882VRM-300E-CH-A</p>
 <p>Suction Trap Assembly (complete)</p>	<p>881VR-300-00-T 881VR-300E-00-T 881VR-300D-00-T 882VR-300-00-T 882VR-300E-00-T 882VR-300D-00-T 882VRM-300-00-T 882VRM-300E-00-T</p>	<p>881VR-300-OH-T 881VR-300E-OH-T 881VR-300D-OH-T 882VR-300-OH-T 882VR-300E-OH-T 882VR-300D-OH-T 882VRM-300-OH-T 882VRM-300E-OH-T</p>	<p>881VR-300-DH-T 881VR-300E-DH-T 881VR-300D-DH-T 882VR-300-DH-T 882VR-300E-DH-T 882VR-300D-DH-T 882VRM-300-DH-T 882VRM-300E-DH-T</p>	<p>881VR-300-DS-T 881VR-300E-DS-T 881VR-300D-DS-T 882VR-300-DS-T 882VR-300E-DS-T 882VR-300D-DS-T 882VRM-300-DS-T 882VRM-300E-DS-T</p>	<p>881VR-300-CH-T 881VR-300E-CH-T 881VR-300D-CH-T 882VR-300-CH-T 882VR-300E-CH-T 882VR-300D-CH-T 882VRM-300-CH-T 882VRM-300E-CH-T</p>
 <p>Trap Fitting Only (No Trap)</p>	<p>881VR-300-00-B 881VR-300E-00-B 881VR-300D-00-B 882VR-300-00-B 882VR-300E-00-B 882VR-300D-00-B 882VRM-300-00-B 882VRM-300E-00-B</p>	<p>881VR-300-OH-B 881VR-300E-OH-B 881VR-300D-OH-B 882VR-300-OH-B 882VR-300E-OH-B 882VR-300D-OH-B 882VRM-300-OH-B 882VRM-300E-OH-B</p>	<p>881VR-300-DH-B 881VR-300E-DH-B 881VR-300D-DH-B 882VR-300-DH-B 882VR-300E-DH-B 882VR-300D-DH-B 882VRM-300-DH-B 882VRM-300E-DH-B</p>	<p>881VR-300-DS-B 881VR-300E-DS-B 881VR-300D-DS-B 882VR-300-DS-B 882VR-300E-DS-B 882VR-300D-DS-B 882VRM-300-DS-B 882VRM-300E-DS-B</p>	<p>881VR-300-CH-B 881VR-300E-CH-B 881VR-300D-CH-B 882VR-300-CH-B 882VR-300E-CH-B 882VR-300D-CH-B 882VRM-300-CH-B 882VRM-300E-CH-B</p>
 <p>DISS Male Vacuum</p>	<p>881VR-300-00-C 881VR-300E-00-C 881VR-300D-00-C 882VR-300-00-C 882VR-300E-00-C 882VR-300D-00-C 882VRM-300-00-C 882VRM-300E-00-C</p>	<p>881VR-300-OH-C 881VR-300E-OH-C 881VR-300D-OH-C 882VR-300-OH-C 882VR-300E-OH-C 882VR-300D-OH-C 882VRM-300-OH-C 882VRM-300E-OH-C</p>	<p>881VR-300-DH-C 881VR-300E-DH-C 881VR-300D-DH-C 882VR-300-DH-C 882VR-300E-DH-C 882VR-300D-DH-C 882VRM-300-DH-C 882VRM-300E-DH-C</p>	<p>881VR-300-DS-C 881VR-300E-DS-C 881VR-300D-DS-C 882VR-300-DS-C 882VR-300E-DS-C 882VR-300D-DS-C 882VRM-300-DS-C 882VRM-300E-DS-C</p>	<p>881VR-300-CH-C 881VR-300E-CH-C 881VR-300D-CH-C 882VR-300-CH-C 882VR-300E-CH-C 882VR-300D-CH-C 882VRM-300-CH-C 882VRM-300E-CH-C</p>
 <p>Tubing Nipple (5/16" Hose Barb)</p>	<p>881VR-300-00-D 881VR-300E-00-D 881VR-300D-00-D 882VR-300-00-D 882VR-300E-00-D 882VR-300D-00-D 882VRM-300-00-D 882VRM-300E-00-D</p>	<p>881VR-300-OH-D 881VR-300E-OH-D 881VR-300D-OH-D 882VR-300-OH-D 882VR-300E-OH-D 882VR-300D-OH-D 882VRM-300-OH-D 882VRM-300E-OH-D</p>	<p>881VR-300-DH-D 881VR-300E-DH-D 881VR-300D-DH-D 882VR-300-DH-D 882VR-300E-DH-D 882VR-300D-DH-D 882VRM-300-DH-D 882VRM-300E-DH-D</p>	<p>881VR-300-DS-D 881VR-300E-DS-D 881VR-300D-DS-D 882VR-300-DS-D 882VR-300E-DS-D 882VR-300D-DS-D 882VRM-300-DS-D 882VRM-300E-DS-D</p>	<p>881VR-300-CH-D 881VR-300E-CH-D 881VR-300D-CH-D 882VR-300-CH-D 882VR-300E-CH-D 882VR-300D-CH-D 882VRM-300-CH-D 882VRM-300E-CH-D</p>

 <p>PB: Puritan-Bennett® 3114-VAC-NT2</p>	 <p>FS 3101-VAC-NT2</p>	 <p>GS 3116-VAC-NT2</p>	 <p>BS 3117-VAC-NT2</p>	 <p>JIS 3108-VAC-NT2</p>	 <p>SIS 3109-VAC-NT2</p>
881VR-300-PB-A 881VR-300E-PB-A 881VR-300D-PB-A 882VR-300-PB-A 882VR-300E-PB-A 882VR-300D-PB-A 882VRM-300-PB-A 882VRM-300E-PB-A	881VR-300-FS-A 881VR-300E-FS-A 881VR-300D-FS-A 882VR-300-FS-A 882VR-300E-FS-A 882VR-300D-FS-A 882VRM-300-FS-A 882VRM-300E-FS-A	881VR-300-GS-A 881VR-300E-GS-A 881VR-300D-GS-A 882VR-300-GS-A 882VR-300E-GS-A 882VR-300D-GS-A 882VRM-300-GS-A 882VRM-300E-GS-A	881VR-300-BS-A 881VR-300E-BS-A 881VR-300D-BS-A 882VR-300-BS-A 882VR-300E-BS-A 882VR-300D-BS-A 882VRM-300-BS-A 882VRM-300E-BS-A	881VR-300-JIS-A 881VR-300E-JIS-A 881VR-300D-JIS-A 882VR-300-JIS-A 882VR-300E-JIS-A 882VR-300D-JIS-A 882VRM-300-JIS-A 882VRM-300E-JIS-A	881VR-300-SIS-A 881VR-300E-SIS-A 881VR-300D-SIS-A 882VR-300-SIS-A 882VR-300E-SIS-A 882VR-300D-SIS-A 882VRM-300-SIS-A 882VRM-300E-SIS-A
881VR-300-PB-T 881VR-300E-PB-T 881VR-300D-PB-T 882VR-300-PB-T 882VR-300E-PB-T 882VR-300D-PB-T 882VRM-300-PB-T 882VRM-300E-PB-T	881VR-300-FS-T 881VR-300E-FS-T 881VR-300D-FS-T 882VR-300-FS-T 882VR-300E-FS-T 882VR-300D-FS-T 882VRM-300-FS-T 882VRM-300E-FS-T	881VR-300-GS-T 881VR-300E-GS-T 881VR-300D-GS-T 882VR-300-GS-T 882VR-300E-GS-T 882VR-300D-GS-T 882VRM-300-GS-T 882VRM-300E-GS-T	881VR-300-BS-T 881VR-300E-BS-T 881VR-300D-BS-T 882VR-300-BS-T 882VR-300E-BS-T 882VR-300D-BS-T 882VRM-300-BS-T 882VRM-300E-BS-T	881VR-300-JIS-T 881VR-300E-JIS-T 881VR-300D-JIS-T 882VR-300-JIS-T 882VR-300E-JIS-T 882VR-300D-JIS-T 882VRM-300-JIS-T 882VRM-300E-JIS-T	881VR-300-SIS-T 881VR-300E-SIS-T 881VR-300D-SIS-T 882VR-300-SIS-T 882VR-300E-SIS-T 882VR-300D-SIS-T 882VRM-300-SIS-T 882VRM-300E-SIS-T
881VR-300-PB-B 881VR-300E-PB-B 881VR-300D-PB-B 882VR-300-PB-B 882VR-300E-PB-B 882VR-300D-PB-B 882VRM-300-PB-B 882VRM-300E-PB-B	881VR-300-FS-B 881VR-300E-FS-B 881VR-300D-FS-B 882VR-300-FS-B 882VR-300E-FS-B 882VR-300D-FS-B 882VRM-300-FS-B 882VRM-300E-FS-B	881VR-300-GS-B 881VR-300E-GS-B 881VR-300D-GS-B 882VR-300-GS-B 882VR-300E-GS-B 882VR-300D-GS-B 882VRM-300-GS-B 882VRM-300E-GS-B	881VR-300-BS-B 881VR-300E-BS-B 881VR-300D-BS-B 882VR-300-BS-B 882VR-300E-BS-B 882VR-300D-BS-B 882VRM-300-BS-B 882VRM-300E-BS-B	881VR-300-JIS-B 881VR-300E-JIS-B 881VR-300D-JIS-B 882VR-300-JIS-B 882VR-300E-JIS-B 882VR-300D-JIS-B 882VRM-300-JIS-B 882VRM-300E-JIS-B	881VR-300-SIS-B 881VR-300E-SIS-B 881VR-300D-SIS-B 882VR-300-SIS-B 882VR-300E-SIS-B 882VR-300D-SIS-B 882VRM-300-SIS-B 882VRM-300E-SIS-B
881VR-300-PB-C 881VR-300E-PB-C 881VR-300D-PB-C 882VR-300-PB-C 882VR-300E-PB-C 882VR-300D-PB-C 882VRM-300-PB-C 882VRM-300E-PB-C	881VR-300-FS-C 881VR-300E-FS-C 881VR-300D-FS-C 882VR-300-FS-C 882VR-300E-FS-C 882VR-300D-FS-C 882VRM-300-FS-C 882VRM-300E-FS-C	881VR-300-GS-C 881VR-300E-GS-C 881VR-300D-GS-C 882VR-300-GS-C 882VR-300E-GS-C 882VR-300D-GS-C 882VRM-300-GS-C 882VRM-300E-GS-C	881VR-300-BS-C 881VR-300E-BS-C 881VR-300D-BS-C 882VR-300-BS-C 882VR-300E-BS-C 882VR-300D-BS-C 882VRM-300-BS-C 882VRM-300E-BS-C	881VR-300-JIS-C 881VR-300E-JIS-C 881VR-300D-JIS-C 882VR-300-JIS-C 882VR-300E-JIS-C 882VR-300D-JIS-C 882VRM-300-JIS-C 882VRM-300E-JIS-C	881VR-300-SIS-C 881VR-300E-SIS-C 881VR-300D-SIS-C 882VR-300-SIS-C 882VR-300E-SIS-C 882VR-300D-SIS-C 882VRM-300-SIS-C 882VRM-300E-SIS-C
881VR-300-PB-D 881VR-300E-PB-D 881VR-300D-PB-D 882VR-300-PB-D 882VR-300E-PB-D 882VR-300D-PB-D 882VRM-300-PB-D 882VRM-300E-PB-D	881VR-300-FS-D 881VR-300E-FS-D 881VR-300D-FS-D 882VR-300-FS-D 882VR-300E-FS-D 882VR-300D-FS-D 882VRM-300-FS-D 882VRM-300E-FS-D	881VR-300-GS-D 881VR-300E-GS-D 881VR-300D-GS-D 882VR-300-GS-D 882VR-300E-GS-D 882VR-300D-GS-D 882VRM-300-GS-D 882VRM-300E-GS-D	881VR-300-BS-D 881VR-300E-BS-D 881VR-300D-BS-D 882VR-300-BS-D 882VR-300E-BS-D 882VR-300D-BS-D 882VRM-300-BS-D 882VRM-300E-BS-D	881VR-300-JIS-D 881VR-300E-JIS-D 881VR-300D-JIS-D 882VR-300-JIS-D 882VR-300E-JIS-D 882VR-300D-JIS-D 882VRM-300-JIS-D 882VRM-300E-JIS-D	881VR-300-SIS-D 881VR-300E-SIS-D 881VR-300D-SIS-D 882VR-300-SIS-D 882VR-300E-SIS-D 882VR-300D-SIS-D 882VRM-300-SIS-D 882VRM-300E-SIS-D

see
previous
page for
complete
table

SUCTION THERAPY

Suction Regulators

Solutions for Life

88XVR Series Suction Regulators (Range 760)

 <p>Vacuum Connector</p> <p>Sucton Connector</p>	 <p>No Adapter 1/8" NPT (F)</p>	 <p>OH: Omeda® 3111-VAC-NT2</p>	 <p>DH: DISS Hand-tight 3222-VAC-NT2</p>	 <p>DS: DISS Hex 3102-VAC-NT2</p>	 <p>CH: Chemetron® 3103-VAC-NT2</p>
 <p>No Adapter 1/8" NPT (F)</p>	<p>882VR-760-00-A 882VR-760E-00-A 882VR-760D-00-A 882VRM-760-00-A 882VRM-760E-00-A</p>	<p>882VR-760-OH-A 882VR-760E-OH-A 882VR-760D-OH-A 882VRM-760-OH-A 882VRM-760E-OH-A</p>	<p>882VR-760-DH-A 882VR-760E-DH-A 882VR-760D-DH-A 882VRM-760-DH-A 882VRM-760E-DH-A</p>	<p>882VR-760-DS-A 882VR-760E-DS-A 882VR-760D-DS-A 882VRM-760-DS-A 882VRM-760E-DS-A</p>	<p>882VR-760-CH-A 882VR-760E-CH-A 882VR-760D-CH-A 882VRM-760-CH-A 882VRM-760E-CH-A</p>
 <p>Suction Trap Assembly (complete)</p>	<p>882VR-760-00-T 882VR-760E-00-T 882VR-760D-00-T 882VRM-760-00-T 882VRM-760E-00-T</p>	<p>882VR-760-OH-T 882VR-760E-OH-T 882VR-760D-OH-T 882VRM-760-OH-T 882VRM-760E-OH-T</p>	<p>882VR-760-DH-T 882VR-760E-DH-T 882VR-760D-DH-T 882VRM-760-DH-T 882VRM-760E-DH-T</p>	<p>882VR-760-DS-T 882VR-760E-DS-T 882VR-760D-DS-T 882VRM-760-DS-T 882VRM-760E-DS-T</p>	<p>882VR-760-CH-T 882VR-760E-CH-T 882VR-760D-CH-T 882VRM-760-CH-T 882VRM-760E-CH-T</p>
 <p>Trap Fitting Only (No Trap)</p>	<p>882VR-760-00-B 882VR-760E-00-B 882VR-760D-00-B 882VRM-760-00-B 882VRM-760E-00-B</p>	<p>882VR-760-OH-B 882VR-760E-OH-B 882VR-760D-OH-B 882VRM-760-OH-B 882VRM-760E-OH-B</p>	<p>882VR-760-DH-B 882VR-760E-DH-B 882VR-760D-DH-B 882VRM-760-DH-B 882VRM-760E-DH-B</p>	<p>882VR-760-DS-B 882VR-760E-DS-B 882VR-760D-DS-B 882VRM-760-DS-B 882VRM-760E-DS-B</p>	<p>882VR-760-CH-B 882VR-760E-CH-B 882VR-760D-CH-B 882VRM-760-CH-B 882VRM-760E-CH-B</p>
 <p>DISS Male Vacuum</p>	<p>882VR-760-00-C 882VR-760E-00-C 882VR-760D-00-C 882VRM-760-00-C 882VRM-760E-00-C</p>	<p>882VR-760-OH-C 882VR-760E-OH-C 882VR-760D-OH-C 882VRM-760-OH-C 882VRM-760E-OH-C</p>	<p>882VR-760-DH-C 882VR-760E-DH-C 882VR-760D-DH-C 882VRM-760-DH-C 882VRM-760E-DH-C</p>	<p>882VR-760-DS-C 882VR-760E-DS-C 882VR-760D-DS-C 882VRM-760-DS-C 882VRM-760E-DS-C</p>	<p>882VR-760-CH-C 882VR-760E-CH-C 882VR-760D-CH-C 882VRM-760-CH-C 882VRM-760E-CH-C</p>
 <p>Tubing Nipple (5/16" Hose Barb)</p>	<p>882VR-760-00-D 882VR-760E-00-D 882VR-760D-00-D 882VRM-760-00-D 882VRM-760E-00-D</p>	<p>882VR-760-OH-D 882VR-760E-OH-D 882VR-760D-OH-D 882VRM-760-OH-D 882VRM-760E-OH-D</p>	<p>882VR-760-DH-D 882VR-760E-DH-D 882VR-760D-DH-D 882VRM-760-DH-D 882VRM-760E-DH-D</p>	<p>882VR-760-DS-D 882VR-760E-DS-D 882VR-760D-DS-D 882VRM-760-DS-D 882VRM-760E-DS-D</p>	<p>882VR-760-CH-D 882VR-760E-CH-D 882VR-760D-CH-D 882VRM-760-CH-D 882VRM-760E-CH-D</p>

 <p>PB: Puritan-Bennett® 3114-VAC-NT2</p>	 <p>FS 3101-VAC-NT2</p>	 <p>GS 3116-VAC-NT2</p>	 <p>BS 3117-VAC-NT2</p>	 <p>JIS 3108-VAC-NT2</p>	 <p>SIS 3109-VAC-NT2</p>
882VR-760-PB-A 882VR-760E-PB-A 882VR-760D-PB-A 882VRM-760-PB-A 882VRM-760E-PB-A	882VR-760-FS-A 882VR-760E-FS-A 882VR-760D-FS-A 882VRM-760-FS-A 882VRM-760E-FS-A	882VR-760-GS-A 882VR-760E-GS-A 882VR-760D-GS-A 882VRM-760-GS-A 882VRM-760E-GS-A	882VR-760-BS-A 882VR-760E-BS-A 882VR-760D-BS-A 882VRM-760-BS-A 882VRM-760E-BS-A	882VR-760-JIS-A 882VR-760E-JIS-A 882VR-760D-JIS-A 882VRM-760-JIS-A 882VRM-760E-JIS-A	882VR-760-SIS-A 882VR-760E-SIS-A 882VR-760D-SIS-A 882VRM-760-SIS-A 882VRM-760E-SIS-A
882VR-760-PB-T 882VR-760E-PB-T 882VR-760D-PB-T 882VRM-760-PB-T 882VRM-760E-PB-T	882VR-760-FS-T 882VR-760E-FS-T 882VR-760D-FS-T 882VRM-760-FS-T 882VRM-760E-FS-T	882VR-760-GS-T 882VR-760E-GS-T 882VR-760D-GS-T 882VRM-760-GS-T 882VRM-760E-GS-T	882VR-760-BS-T 882VR-760E-BS-T 882VR-760D-BS-T 882VRM-760-BS-T 882VRM-760E-BS-T	882VR-760-JIS-T 882VR-760E-JIS-T 882VR-760D-JIS-T 882VRM-760-JIS-T 882VRM-760E-JIS-T	882VR-760-SIS-T 882VR-760E-SIS-T 882VR-760D-SIS-T 882VRM-760-SIS-T 882VRM-760E-SIS-T
882VR-760-PB-B 882VR-760E-PB-B 882VR-760D-PB-B 882VRM-760-PB-B 882VRM-760E-PB-B	882VR-760-FS-B 882VR-760E-FS-B 882VR-760D-FS-B 882VRM-760-FS-B 882VRM-760E-FS-B	882VR-760-GS-B 882VR-760E-GS-B 882VR-760D-GS-B 882VRM-760-GS-B 882VRM-760E-GS-B	882VR-760-BS-B 882VR-760E-BS-B 882VR-760D-BS-B 882VRM-760-BS-B 882VRM-760E-BS-B	882VR-760-JIS-B 882VR-760E-JIS-B 882VR-760D-JIS-B 882VRM-760-JIS-B 882VRM-760E-JIS-B	882VR-760-SIS-B 882VR-760E-SIS-B 882VR-760D-SIS-B 882VRM-760-SIS-B 882VRM-760E-SIS-B
882VR-760-PB-C 882VR-760E-PB-C 882VR-760D-PB-C 882VRM-760-PB-C 882VRM-760E-PB-C	882VR-760-FS-C 882VR-760E-FS-C 882VR-760D-FS-C 882VRM-760-FS-C 882VRM-760E-FS-C	882VR-760-GS-C 882VR-760E-GS-C 882VR-760D-GS-C 882VRM-760-GS-C 882VRM-760E-GS-C	882VR-760-BS-C 882VR-760E-BS-C 882VR-760D-BS-C 882VRM-760-BS-C 882VRM-760E-BS-C	882VR-760-JIS-C 882VR-760E-JIS-C 882VR-760D-JIS-C 882VRM-760-JIS-C 882VRM-760E-JIS-C	882VR-760-SIS-C 882VR-760E-SIS-C 882VR-760D-SIS-C 882VRM-760-SIS-C 882VRM-760E-SIS-C
882VR-760-PB-D 882VR-760E-PB-D 882VR-760D-PB-D 882VRM-760-PB-D 882VRM-760E-PB-D	882VR-760-FS-D 882VR-760E-FS-D 882VR-760D-FS-D 882VRM-760-FS-D 882VRM-760E-FS-D	882VR-760-GS-D 882VR-760E-GS-D 882VR-760D-GS-D 882VRM-760-GS-D 882VRM-760E-GS-D	882VR-760-BS-D 882VR-760E-BS-D 882VR-760D-BS-D 882VRM-760-BS-D 882VRM-760E-BS-D	882VR-760-JIS-D 882VR-760E-JIS-D 882VR-760D-JIS-D 882VRM-760-JIS-D 882VRM-760E-JIS-D	882VR-760-SIS-D 882VR-760E-SIS-D 882VR-760D-SIS-D 882VRM-760-SIS-D 882VRM-760E-SIS-D

see
previous
page for
complete
table

SUCTION THERAPY

Suction Regulator

Solutions for Life



880VS Series Venturi Suction Regulator

Glow in the Dark Gauge

- 2" gauge, easy to read
- Color coded flow ranges
- Dual Scale mmHg and kPa

Safety Device

- Positive pressure relief valve provides safe use
- ISO 10079-3 Compliant

Control Knob

- Easy adjustment of vacuum pressure
- No Leakage

Suction Trap

- Prevents contamination caused by overflow



Type of Adapter

- None: Without adapter
- OH: Ohmeda® adapter
- DS: DISS hex adapter
- DH: DISS Handwheel adapter
- CH: Chemetron® adapter
- PB: Puritan-Bennett® adapter
- FS: French (NF S 90-116) adapter
- GS: German (DIN 13260-2) adapter
- BS: British (BS5682-1998) adapter
- JIS: Japanese Style adapter
- SIS: Australian (AS2896) adapter

Venturi Suction

- Delivers steady and high vacuum flow from a medical air or oxygen gas source
- Control Range: -1 kPa ~ -80 kPa

Suction Connector

- GENTEC® provides a wide range of fittings to meet customers' needs

880VS-760-00-A-T

Technical Specification

Dimension (H x W x D)	90 x 90 x 60 mm
Weight	490 g
Range	-1 kPa ~ -80 kPa
Maximum Suction	> 25 L / min
Source of Pressure	300 ~ 500 kPa
Vacuum Control	Needle valve control venturi
Vacuum Gauge Range	0 ~ -100 kPa / 0 ~ 760 mmHg

Ordering Information

Ex: 880VS	- 760	- 00	- A	- T
Series	Measuring Range	Vacuum Connection	Driving Gas	Suction Connection
880VS	760: 0~760 mmHg	00: No adapter, 1/8" NPT (F) OH: 3111U (Ohmeda) adapter DH: 3222U (DISS hand-tight) adapter DS: 3102 (DISS Hex) adapter CH: 3103 (Chemetron) adapter PB: 3104 (Puritan-Bennett) adapter FS: 3105 (FS) adapter GS: 3116 (GS) adapter BS: 3117 (BS) adapter JIS: 3108 (JIS) adapter SIS: 3109 (SIS) adapter	A: Medical Air O: Oxygen	00: No adapter, 1/8" NPT (F) T: With suction trap assembly (complete) TD: With suction trap (complete DISS Female nut) TH: With suction trap (complete DISS female vacuum hand-tight) H: With tubing nipple (5/16" hose barb) D: With DISS male vacuum N: With NIST adapter F: With trap fitting only (no trap)

SUCTION THERAPY

Suction Regulator

Solutions for Life

880VS Series Venturi Suction Regulator

 Vacuum Connector Suction Connector	 No Adapter No Adapter, 1/8" NPT (F)	 OH: Omeda® 3101-OXY-NT2 3111-AIR-NT2	 DH: DISS Hand-tight 3222-OXY-NT2 3222-AIR-NT2	 DS: DISS Hex 3102-OXY-NT2 3102-AIR-NT2	 CH: Chemetron® 3103-OXY-NT2 3103-AIR-NT2
 No Adapter No Adapter, 1/8" NPT (F)	880VS-760-00-A-00 880VS-760-00-O-00	880VS-760-OH-A-00 880VS-760-OH-O-00	880VS-760-DH-A-00 880VS-760-DH-O-00	880VS-760-DS-A-00 880VS-760-DS-O-00	880VS-760-CH-A-00 880VS-760-CH-O-00
 Suction Trap Assembly(complete)	880VS-760-00-A-T 880VS-760-00-O-T	880VS-760-OH-A-T 880VS-760-OH-O-T	880VS-760-DH-A-T 880VS-760-DH-O-T	880VS-760-DS-A-T 880VS-760-DS-O-T	880VS-760-CH-A-T 880VS-760-CH-O-T
 Suction Trap (with DISS female nut & nipple vacuum fitting)	880VS-760-00-A-TD 880VS-760-00-O-TD	880VS-760-OH-A-TD 880VS-760-OH-O-TD	880VS-760-DH-A-TD 880VS-760-DH-O-TD	880VS-760-DS-A-TD 880VS-760-DS-O-TD	880VS-760-CH-A-TD 880VS-760-CH-O-TD
 Suction Trap (with DISS female vacuum hand-tight)	880VS-760-00-A-TH 880VS-760-00-O-TH	880VS-760-OH-A-TH 880VS-760-OH-O-TH	880VS-760-DH-A-TH 880VS-760-DH-O-TH	880VS-760-DS-A-TH 880VS-760-DS-O-TH	880VS-760-CH-A-TH 880VS-760-CH-O-TH
 Tubing Nipple (5/16" Hose Barb)	880VS-760-00-A-H 880VS-760-00-O-H	880VS-760-OH-A-H 880VS-760-OH-O-H	880VS-760-DH-A-H 880VS-760-DH-O-H	880VS-760-DS-A-H 880VS-760-DS-O-H	880VS-760-CH-A-H 880VS-760-CH-O-H
 DISS Male Vacuum	880VS-760-00-A-D 880VS-760-00-O-D	880VS-760-OH-A-D 880VS-760-OH-O-D	880VS-760-DH-A-D 880VS-760-DH-O-D	880VS-760-DS-A-D 880VS-760-DS-O-D	880VS-760-CH-A-D 880VS-760-CH-O-D
 NIST vacuum fitting	880VS-760-00-A-N 880VS-760-00-O-N	880VS-760-OH-A-N 880VS-760-OH-O-N	880VS-760-DH-A-N 880VS-760-DH-O-N	880VS-760-DS-A-N 880VS-760-DS-O-N	880VS-760-CH-A-N 880VS-760-CH-O-N
 Trap Fitting Only (No Trap)	880VS-760-00-A-F 880VS-760-00-O-F	880VS-760-OH-A-F 880VS-760-OH-O-F	880VS-760-DH-A-F 880VS-760-DH-O-F	880VS-760-DS-A-F 880VS-760-DS-O-F	880VS-760-CH-A-F 880VS-760-CH-O-F

 <p>PB: Puritan-Bennett® 3104-OXY-NT2 3114-AIR-NT2</p>	 <p>FS 3101-OXY-NT2 3101-AIR-NT2</p>	 <p>GS 3116-OXY-NT2 3116-AIR-NT2</p>	 <p>BS 3117-OXY-NT2 3117-AIR-NT2</p>	 <p>JIS 3108-OXY-NT2 3108-AIR-NT2</p>	 <p>SIS 3109-OXY-NT2 3109-AIR-NT2</p>
880VS-760-PB-A-00 880VS-760-PB-O-00	880VS-760-FS-A-00 880VS-760-FS-O-00	880VS-760-GS-A-00 880VS-760-GS-O-00	880VS-760-BS-A-00 880VS-760-BS-O-00	880VS-760-JIS-A-00 880VS-760-JIS-O-00	880VS-760-SIS-A-00 880VS-760-SIS-O-00
880VS-760-PB-A-T 880VS-760-PB-O-T	880VS-760-FS-A-T 880VS-760-FS-O-T	880VS-760-GS-A-T 880VS-760-GS-O-T	880VS-760-BS-A-T 880VS-760-BS-O-T	880VS-760-JIS-A-T 880VS-760-JIS-O-T	880VS-760-SIS-A-T 880VS-760-SIS-O-T
880VS-760-PB-A-TD 880VS-760-PB-O-TD	880VS-760-FS-A-TD 880VS-760-FS-O-TD	880VS-760-GS-A-TD 880VS-760-GS-O-TD	880VS-760-BS-A-TD 880VS-760-BS-O-TD	880VS-760-JIS-A-TD 880VS-760-JIS-O-TD	880VS-760-SIS-A-TD 880VS-760-SIS-O-TD
880VS-760-PB-A-TH 880VS-760-PB-O-TH	880VS-760-FS-A-TH 880VS-760-FS-O-TH	880VS-760-GS-A-TH 880VS-760-GS-O-TH	880VS-760-BS-A-TH 880VS-760-BS-O-TH	880VS-760-JIS-A-TH 880VS-760-JIS-O-TH	880VS-760-SIS-A-TH 880VS-760-SIS-O-TH
880VS-760-PB-A-H 880VS-760-PB-O-H	880VS-760-FS-A-H 880VS-760-FS-O-H	880VS-760-GS-A-H 880VS-760-GS-O-H	880VS-760-BS-A-H 880VS-760-BS-O-H	880VS-760-JIS-A-H 880VS-760-JIS-O-H	880VS-760-SIS-A-H 880VS-760-SIS-O-H
880VS-760-PB-A-D 880VS-760-PB-O-D	880VS-760-FS-A-D 880VS-760-FS-O-D	880VS-760-GS-A-D 880VS-760-GS-O-D	880VS-760-BS-A-D 880VS-760-BS-O-D	880VS-760-JIS-A-D 880VS-760-JIS-O-D	880VS-760-SIS-A-D 880VS-760-SIS-O-D
880VS-760-PB-A-N 880VS-760-PB-O-N	880VS-760-FS-A-N 880VS-760-FS-O-N	880VS-760-GS-A-N 880VS-760-GS-O-N	880VS-760-BS-A-N 880VS-760-BS-O-N	880VS-760-JIS-A-N 880VS-760-JIS-O-N	880VS-760-SIS-A-N 880VS-760-SIS-O-N
880VS-760-PB-A-F 880VS-760-PB-O-F	880VS-760-FS-A-F 880VS-760-FS-O-F	880VS-760-GS-A-F 880VS-760-GS-O-F	880VS-760-BS-A-F 880VS-760-BS-O-F	880VS-760-JIS-A-F 880VS-760-JIS-O-F	880VS-760-SIS-A-F 880VS-760-SIS-O-F

see
previous
page for
complete
table

SUCTION THERAPY

Accessories

Solutions for Life

Suction Traps

Suction traps prevent contamination caused by overflow

880VT-0	880VT-1	880VT-2	880VT-3
			
with no fitting on top (1/8" NPT female)	with male trap fitting only	with DISS female nut & nipple vacuum fitting	with DISS female vacuum hand-tight

880VT-4	880VT-5	880VT-6	880VT-7	880VT-8
				
with tubing nipple (5/16" hose barb)	with complete trap fitting assembly (male & female)	with complete DISS Female nut & nipple vacuum fitting	with complete DISS female vacuum hand-tight	with small filter, disposable

Suction Gauges

<p>Clockwise</p>	 <p>GR25PM-V160 (0 ~ 160 mmHg)</p>	 <p>GR25PM-V200 (0 ~ 200 mmHg)</p>	 <p>GR25PM-V300 (0 ~ 300 mmHg)</p>	 <p>GR25PM-V760 (0 ~ 760 mmHg)</p>
<p>Counter-Clockwise</p>	 <p>GR25PM-V160-E (0 ~ 160 mmHg)</p>	 <p>GR25PM-V300-E (0 ~ 300 mmHg)</p>	 <p>GR25PM-V760-E (0 ~ 760 mmHg)</p>	

Venturi Vacuum Modules

Features

- Designed to work with vacuum regulator
- Small and compact
- With positive pressure relief valve
- With silencer
- Maximum vacuum: 0 ~ 600 mmHg
- Maximum air free flow: 30 L / min
- Noise: < 60 dBa



880VE-X




SUCTION THERAPY

Medical Ejector Suction System

Solutions for Life

Venturi System



Ordering Information

Ex: 881VR - 300	E	- VG	A	- 00	- T	- 1L	* 2	
Series	Measuring Range	Pointer Rotate Mode	Suction Generator	Driving Gas	Vacuum Connection	Suction Connection	Canister Volumn	Canister Quantity
881VR	160: 0~160 mmHg	None:	VG: With Venturi	A: Air	00: No adapter	T: Suction trap assembly (complete)	1L: 1L canister	1: one canister
882VR	200: 0~200 mmHg 300: 0~300 mmHg 760: 0~760 mmHg (only for 882VR)	Clockwise, (USA) E: Counter-clockwise (ISO)	Vacuum module None: Without Venturi Vacuum module	O: Oxygen	1/8" NPT (F) FDS: DISS hose barb (USA) FNS: NIST hose barb (ISO)	None: Without suction trap	1.5L: 1.5L canister 2.5L: 2.5L canister	2: two canisters

Standard suction hose length is 0.5 m, for other length, please contact GENTEC.


Disposable Collection

- Constructed with durable polyethylene material
- Used with polycarbonate reusable collection jars
- Liner can be easily disposed and replaced
- Tandem port for multiple canister connections
- Overflow protection and non-return valves prevent backflow
- Can be used with **GENTEC** re-usable collection jar
- Holds 2 liters
- Hose barb inlet and outlet



NEW

Use of Liners




SUCTION THERAPY

Accessories


Solutions for Life

Disposable Collection

Two Liners Connected in Tandem



Four Liners Connected in Tandem



Four Liners Connected in Tandem

Ordering Information

DL-2000B	DL-2000W	615-901	615-902	615-903	615-904	615-905
Liner with overflow protection, 2L (Φ 6-8 hose barb)	Liner without overflow protection, 2L (Φ 6-8 hose barb)	Collections Jar, 2L, reusable	Suction hose, 0.35 m with suction tip	Suction hose, 2 m with suction tip	Suction hose, 3 m with suction tip	Suction hose, 10 m without suction tip


Reusable Collection

SJ Series Reusable Jars

- 1 liter, 1.5 liter, and 2.5 liter models available
- The SJ-1000 series comes with a polycarbonate jar and anodized aluminum cover
- The SJ-1500 and SJ-2500 Series comes with polycarbonate jar and cover
- Built-in bracket for use with GENTEC® wall slides
- Polycarbonate jar for maximum durability
- Ball float design prevents overflow
- GENTEC® SJ series reusable canisters are manufactured to be durable and autoclavable (up to 121°C)

DJ Series Disposable Jars

- 0.85 liter, 1 liter, and 2 liter models available
- comes with polycarbonate jar and cover
- Built-in bracket for use with GENTEC® wall slides
- Ball float design prevents overflow
- Hose barb inlet and outlet



Ordering Information

Model	SJ Series				DJ Series		
							
Volume	1 L	1 L	1.5 L	2.5 L	0.85 L	1 L	2 L
Overfill Protection Devices	✓	✓	✓	✓	✓	✓	✓
Vacuum Connection	DISS	Hose Barb (φ 6-8)	Hose Barb (φ 6-8)	Hose Barb (φ 6-8)	Hose Barb (φ 6-8)	Hose Barb (φ 6-8)	Hose Barb (φ 6-8)
Patient Connection	Hose Barb (φ 6-8)						

SUCTION THERAPY

Accessories

Solutions for Life

Wall Slides

- For use on all surfaces
- Fits GENTEC DJ, SJ & DL series canister



SJ-CLP



CLP-D40-45



3910-Sld



3820-Sld



Collection Stations

- Designed to support collection jars
- Stable 5-caster base with two locking casters
- Durable construction, easy to clean
- Supports up to 4 canisters
- Includes suction trap for additional safety
- Shut-off valve and canister slides included



810-00(DH) / 810-00(HB)

Suction Tubings

Model	Description
ST6-4	Suction Tubing 1/4" (6mm), I.D. 4 Meters, Without Suction Tip
ST6-100	Suction Tubing 1/4" (6mm), I.D. 100 Meters, Without Suction Tip
ST8-4	Suction Tubing 5/16" (8mm), I.D. 4 Meters, With Suction Tip
ST8-160	Suction Tubing 5/16" (8mm), I.D. 160 Meters, With Suction Tip



Disposable Bacteria Filters

- Prevents contamination and fluid from entering the vacuum regulator and vacuum systems
- PVDF special filter material, can effectively filter the virus and bacteria
- Effective filtration area: 32 cm²
- Maximum operating temperature: 140 °F (60 °C)
- Maximum operating pressure: -14.5 ~ 14.5 psi
- Filters down to 0.22 micron




F11-HB



F11-NT2/HB



Dimensions



F11-HB

F11-NT2/HB

Ordering Information

Model	Description
F11-NT2/HB-0.45	Hydrophobic, Inlet: 1/4" Hose Barb, Outlet: 1/8" NPT, Filter Accuracy: 0.45µm
F11-NT2/HB-0.22	Hydrophobic, Inlet: 1/4" Hose Barb, Outlet: 1/8" NPT, Filter Accuracy: 0.22µm
F11-HB-0.45	Hydrophobic, Inlet and Outlet: 1/4" Hose Barb, Filter Accuracy: 0.45µm
F11-HB-0.22	Hydrophobic, Inlet and Outlet: 1/4" Hose Barb, Filter Accuracy: 0.22µm